

INSTALLING PERIOPAL ON WINDOWS 7

We like to think of the Windows 7 and the XP Mode Virtual Machine as two distinct computers. PerioPal should NOT be installed on the Windows 7 computer but only in the XP Mode computer. Please follow all steps and their links.



1 Install XP Mode in Windows 7:

- After installing XP Mode, the icon will be seen on the Windows 7 desktop.
- Click the Windows XP icon to open the Windows XP Mode
- Find the **internet browser** in XP Mode to download Adobe Reader and PerioPal

2 In XP Mode, [download](#) and install Adobe Reader. [Configure](#) Adobe Reader for PerioPal.

3 In XP Mode, download the [Free Download Version](#) of PerioPal:.

- Single PC Download
- Network Install [See instructions](#)
- Register your 30-Day Trial of PerioPal
- [Create a Desktop icon for PerioPal](#)

4 Order the USB Buddy Sound Device and Andrea ANC 750 Headset Combination.

- This combination has been tested and is successfully recognized by the Device Manager in Windows 7-XP Mode.

5 Configure and attach the Buddy USB and ANC 750 Microphone

- Return to XP Mode. Attach the USB Buddy in XP Mode.
- **Learn to exit the XP Mode properly** – keeping XP in 'Hibernation'.
- On the Windows 7 Desktop, attach the USB Buddy in the Device Manager.
* * * * *
- Using the PerioPal icon on the Windows 7 desktop, open PerioPal.
- Configure the Audio Settings.
- Configure the Speech Settings.
- Create and train a Voice Profile.
- Set your User Preferences
- Enter a Patient



6 Learn how to launch PerioPal everyday:

- Print and laminate: [Routine Maintenance for Windows 7](#).
- Make a copy for each workstation.
- Understand that Windows 7 does not 'remember' select screens when closed.
- Understand you will not have proper sound input unless select screens are correct.

7 Your printer is installed on Windows 7. Reports are printed in Windows 7 – not in XP Mode.

8 [Network PerioPal](#). Contact Support for Bridge options.

- [Create a Shared Folder and Map a Network Drive](#)

UNINSTALL PERIOPAL ON WINDOWS 7

1 PerioPal should not be installed on the Windows 7 computer but only on the XP Mode computer.

[Uninstalling PerioPal from the XP machine is identical to steps used in Windows XP .](#)

- Control Panel | Add Remove Programs – uninstall SAPI Voice Tools and PPpdf.(nova)
- Delete the folder C:\PPserver
- Delete the folder C:\Program Files\PerioPal
- **If PerioPal is networked**, delete the mapped drive to the server folder and the server folder C:\PPserver.

2 If PerioPal is installed on Windows 7 the [process to remove](#) is exactly the same, except that the Control Panel application has been renamed “**Programs and Features**”

UNINSTALL XP MODE ON WINDOWS 7

1 Windows 7 | Control Panel | Programs and Features – uninstall Windows XP Mode.