

Download, Install, and Register the Single PC Version:

1 On the Home page, click **DOWNLOAD NOW** button.

2 On the Download page, click the **Single PC Install** button to begin the download.



- You have reviewed the Hardware Requirements.
- You may have to temporarily disable your anti-virus software.
- You have completely [removed previous installations](#) of PerioPal.


3 The file will download. When complete, click the **setup.exe** file to start the installation.




4 You will be asked for permission to install. PerioPal is a secure website. Click **OK**.




5 The Install Wizard Screens prepare to install. Please be patient.




6 The Installation Wizard Screens begin. Follow the directions by clicking Next.




7 The Workstation Installation is complete. Click **Finish**, but do not close. The NOVA PRINTER and SAPI VOICE TOOLS will install.



8 Return to your desktop. Click the PerioPal icon to register your software.




9 The **Registration** Screen will appear. Enter your practice information as you would like it to appear on Reports. Click **Register On-line**. View Privacy Policy.



10 The Training Screen will appear. Close this screen to find the Preferences Screen. You can also find the **Training Materials** link on PerioPal's Download page.



11 The Preferences Screen:



- **Registered As:** Your System ID and Registration information.
- **Subscription Box** displays you Approval Date and Order Online option.

Troubleshooting Links:

[The Program Downloads But Will Not Register](#)

Error With Installation or Registration?

[Remove PerioPal](#) and try again.

[With my second installation, PerioPal will not let me register.](#)